## ngs\_week3 assignment - Xiaoyan Wen

XiaoyanWen

2/6/2022

## Next Generation Sequence Analysis Homework Week 3

In this assignment, you will align short reads to the human reference genome and conduct a short exercise that uses samtools to process SAM/BAM short read alignments. This is part of multi-week "re-sequencing" workflow where the aim is to call single nucleotide polymorphisms (SNPs) from Illumina short read sequencing of human genomes.

##About the data This week we continue working with re-sequencing data from the 1000 Genomes Project. We will work with 30 adapter-trimmed and quality filtered reads (processed via fastp as in Week 2).

##Task 1: Preparing the human reference genome ###Introduction to reference genomes Reference genomes are typically FASTA-formatted genome sequences representing the primary sequence of a (typically) haploid copy of the genome. In this task, you will prepare a copy human reference genome for short read alignment with the Burrows-Wheeler Aligner BWA-MEM.

Choosing the appropriate human reference genome is complicated by the fact that many versions exist (GRCh38/hg38 is current) and for each version hundreds of variants exists. For example, there are repeat-masked, soft-masked, primary assembly and top-level assemblies (with some variable regions appearing more than once as different haplotypes, or alternate sequences) to name just a few.

we chose an appropriate variant of the hg38 version human reference genome for short read alignment from ENSEMBL. The FASTA-formatted file was then downloaded using

```
wget ftp://ftp.ensembl.org/pub/current_fasta/homo_sapiens/dna/Homo_sapiens.GRCh38.dna_sm.primary_assembly.fa.gz
```

###Normalizing the sequence identifiers of a reference genome One disadvantage of FASTA format is that it does not include strict rules for sequence identifiers. Many downstream tools require that the sequence identifiers contain no white spaces (even though the FASTA specification allows whitespaces in identifiers).

####Picard-Tools NormalizeFasta program: https://broadinstitute.github.io/picard/command-line-overview.htm/#NormalizeFasta (https://broadinstitute.github.io/picard/command-line-overview.htm/#NormalizeFasta) This tool has been used to removes any white spaces in sequence identifiers in the reference FASTA to avoid possible conflicts with downstream applications.

###Index files for the reference genome FASTA Many NGS software require that input files be indexed prior to including them in work flows. Index files permit rapid lookups of genome coordinates in large files and dramatically speed up computation times.

The re-sequencing workflow we will run to align reads from samples in the 1,000 genomes project followed by snp-calling requires the following index files: - a FASTA index - a set of index files required by the BWA-MEM aligner - a dictionary file required by the Genome Analysis Toolkit (GATK)

Below, you will create a FASTA index and a bwa index for the normalized reference genome. We will create a the dictionary file in a subsequence exercise.

-See here for requirements for the reference by GATK: https://gatk.broadinstitute.org/hc/en-us/articles/360035531652-FASTA-Reference-genome-format (https://gatk.broadinstitute.org/hc/en-us/articles/360035531652-FASTA-Reference-genome-format) -See samtools manual for the samtools faidx command: http://www.htslib.org/doc/samtools.html) -See BWA manual for bwa index here: http://bio-bwa.sourceforge.net/bwa.shtml (http://bio-bwa.sourceforge.net/bwa.shtml)

Before beginning, log in to Greene and request a compute node for interactive use:

```
srun --time=4:00:00 --mem=4GB --pty /bin/bash
```

Create a directory in your /scratch as follows, change directories to it, and copy to it the normalized hg38 reference genome FASTA as follows:

```
cd $SCRATCH
mkdir hg38
cd hg38
cp /scratch/work/courses/BI7653/hw3.2022/hg38/Homo_sapiens.GRCh38.dna_sm.primary_assembly.normalized.fa .
```

Create a slurm job script that will create the FASTA index files in the same directory as the reference genome FASTA. You may use the template slurm script provided by your instructor in Week 1

```
/scratch/work/courses/BI7653/hw1.2022/slurm_template.sh
```

Key point: For virtually all NGS applications, the reference genome index files must appear in the same directory as the reference FASTA. Your slurm job script should do the following: - Request 32GB of memory for 5 hours. - Load the samtools and bwa modules (most recent versions on Greene) - Run samtools faidx and bwa index using the following commands on the normalized reference fasta you copied to your hg38 directory above.

```
samtools faidx <reference fasta>
bwa index -a bwtsw <reference fasta>
```

You may now execute your job script. This should take less than 2 hours from the time the job is released from the queue. Occasionally, students have encountered "segmentation fault" errors when running this tool. If your STDERR for the job contains an error of this type, try running again.

Q1.1. Please report the contents of your job script [ 1 point ].

```
#!/bin/bash
#
#SBATCH --nodes=1
#SBATCH --tasks-per-node=1
#SBATCH --true=5:00:00
#SBATCH --true=5:00:00
#SBATCH --true=5:00:00
#SBATCH --mem=32GB
#SBATCH --mem=32GB
#SBATCH --mail-type=FAIL
#SBATCH --mail-user=xw2470@nyu.edu
module purge
echo script begin: $(date)
module load samtools/intel/1.14
module load bwa/intel/0.7.17
samtools faidx Homo_sapiens.GRCh38.dna_sm.primary_assembly.normalized.fa
bwa index -a bwtsw Homo_sapiens.GRCh38.dna_sm.primary_assembly.normalized.fa
echo script completed: $(date)
```

Q1.2. Upon job completion, please execute Is -al in your hg38 directory and report the output [ 1 point ]. screenshot:

## Task 2: Short read alignment with BWA-MEM

BWA-MEM is a common alignment tool that was used in the 1000 Genomes Project and is the preferred tool in the Genome Analysis Toolkit SNP-calling and genotyping workflow which we will follow in this course. You may wish to read the BWA documentation http://bio-bwa.sourceforge.net/bwa.shtml (http://bio-bwa.sourceforge.net/bwa.shtml)

In this task, you will execute a job script provided by your instructor to align paired-end fastq data from 30 human samples to the human reference genome. Your instructor has prepared a set of 30 samples (60 paired-end fastqs) for you to align (paired end sequence data from the same 21 samples processed with fastp in week 2 + 9 additional paired end samples) representing four populations in the 1000 genomes project. The fastqs were processed using the same fastp command line you used in week 2 Task 2 and are ready for short read alignment with BWA-MEM.

You will now align short reads to the human reference genome using a script provided by your instructor. You are simply asked to modify one variable in the script to contain the path on your /scratch to the reference genome FASTA (see below)

Create a ngs.week3 directory and Task 2 subdirectory in your /scratch directory.

```
cd $SCRATCH
pwd
mkdir ngs.week3
cd ngs.week3
mkdir task2
cd task2
cd task2
cp /scratch/work/courses/BI7653/hw3.2022/hw3_bwamem.slurm . # copy script to present workding directory
```

The script executes a job array that will find the paired-end fastqs to be aligned using a tab-delimited table with columns sample name, read 1 fastq file name, and read 2 fastq filename. You may review that file here (e.g., with less command):

```
/scratch/work/courses/BI7653/hw3.2022/fastqs.processed/hw3_fastqs.processed.txt
```

The processed fastqs are located here:

```
/scratch/work/courses/BI7653/hw3.2022/fastqs.processed
```

We will align the reads with BWA-MEM using the following settings:

```
bwa mem -M -t $$LURM_CPUS_PER_TASK -R "@RG\tID:${sample}.id\tSM:${sample}\tPL:ILLUMINA\tLB:${sample}.lb" <reference fasta> < read1.fq> <read2.fq> > <sam>
```

Note the following: -The SLURM\_CPUS\_PER\_TASK environmental variable will take on the value of 8 because the slurm directives at the top of the script used #SBATCH -cpus-per-task=8. -The -R argument adds an @RG read group header line to the output SAM file. This gives the sample unique identifer and information tags in the output SAM header.

Now review hw3\_bwamem.slurm script which you copied to your /scratch directory above. Convince yourself that #SBATCH array=1-30 directive and subsequent code will correctly process the 30 samples listed in the file

```
/scratch/work/courses/BI7653/hw3.2022/fastqs.processed/hw3_fastqs.processed.txt
```

Now modify the script to specify the ref variable. The ref variable should be set to the full path (from root = "/") to the hg38 reference genome fasta that you indexed in Task 1. If you need to, you may cd to the directory with the reference genome and use pwd to view the path to that directory from root.

Verify that you are in your Task 2 directory and execute the job script:

```
sbatch hw3_bwamem.slurm
```

Execute the following immediately after execution to check on the status of your job:

```
squeue -u xw2470
```

```
Q2.1 Now either take a screen shot showing your squeue command and the output (or copy the output to your homework report) [1 point]

[xx/2470@cs011 task2]$ sbatch hw3_bwamem.slurm

[xw/2470@cs011 task2]$ sequeue -u xw/2470

[xw/2470@cs011 task2]$ sequeue -u xw/2470
  oash: sequeue: command not found
[xw2470@cs011 task2]$ squeue -u xw2470
JOBID PARTITION NAME
                                                                                                                                    TIME
0:37
0:37
0:37
                                                                                                                                                   NODES NODELIST(REASON)
                  14853877_1
14853877_2
14853877_3
14853877_4
                                                                                                                                                            1 cs029
1 cs029
1 cs050
                                                             cs bwamem_a
                                                                                               xw2470
                                                                                                                 cs bwamem_a
cs bwamem_a
cs bwamem_a
cs bwamem_a
                                                                                               xw2470
xw2470
                                                                                                                                    0:37
0:37
0:37
                                                                                               xw2470
                                                                                                                                                             1 cs057
1 cs059
                   14853877_4
14853877_5
14853877_6
14853877_7
14853877_8
14853877_9
                                                                                               xw2470
xw2470
                                                                                                                                                             1 cs059
1 cs061
                                                             cs bwamem_a
cs bwamem_a
cs bwamem_a
                                                                                               xw2470
xw2470
xw2470
                                                                                                                                    0:37
0:37
0:37
                                                                                                                                                             1 cs064
1 cs070
                                                             cs bwamem_a
cs bwamem_a
cs bwamem_a
                                                                                                                                                             1 cs072
                  14853877_10
14853877_11
14853877_12
                                                                                               xw2470
xw2470
                                                                                                                                    0:37
0:37
                                                                                                                                                             1 cs072
1 cs080
                                                            cs bwamem_a
cs bwamem_a
cs bwamem_a
cs bwamem_a
cs bwamem_a
                                                                                              xw2470
xw2470
xw2470
xw2470
xw2470
                                                                                                                                    0:37
0:37
0:37
                  14853877_13
14853877_14
                                                                                                                                                             1 cs082
1 cs082
                  14853877_15
14853877_16
14853877_17
                                                                                                                                    0:37
0:37
0:37
                                                             cs bwamem_a
                                                                                               xw2470
                 14853877_18
14853877_19
14853877_20
14853877_21
14853877_21
                                                                                               xw2470
xw2470
xw2470
                                                                                                                                    0:37
0:37
0:37
                                                             cs bwamem_a
cs bwamem_a
                                                                                                                                                             1 cs099
1 cs106
                                                             cs bwamem_a
cs bwamem_a
cs bwamem_a
                                                                                                                                                             1 cs122
                                                                                               xw2470
xw2470
                                                                                                                                                             1 cs125
1 cs130
                  14853877_23
14853877_24
14853877_25
                                                                                                                                    0:37
0:37
0:37
                                                                                                                                                             1 cs133
1 cs146
1 cs146
                                                             cs bwamem_a
cs bwamem_a
                                                                                               xw2470
                                                                                               xw2470
xw2470
                                                             cs bwamem_a
                                                                                                                                    0:37
0:37
0:37
                                                                                                                                                             1 cs146
1 cs149
1 cs149
                  14853877_26
14853877_27
                                                                                               xw2470
xw2470
                                                             cs bwamem_a
cs bwamem_a
                  14853877 28
                                                                                               xw2470
                                                                                                                                    0:37
0:37
2:51
                  14853877_29
14853877_30
                                                             cs bw
                                                                                               xw2470
xw2470
                                                                                                                                                             1 cs157
1 cs160
                                                                              bash
                                                                                               xw2470
                                                                                                                                                             1 cs011
                        011 task2]$
```

This job could take 4-8 hours to complete once each of the 30 "subjobs" have executed. When your job is complete, please return to the directory where you executed the script and report the output of the following command.

```
grep _ESTATUS_ slurm-14853877*.out
```

```
w2470@gr062 task2]$ grep _ESTATUS_ slurm-14853877*.o
                                                                      bwa mem for HG00149
bwa mem for HG00260
                                                                      bwa mem for NA18907
bwa mem for NA19137
                                                                     wa mem for NA19137
bwa mem for NA19093
bwa mem for NA19256
bwa mem for NA19098
bwa mem foe
                                                                      bwa mem for NA18870
bwa mem for NA18909
                                                                    bwa mem for NA19138
bwa mem for NA18757 ]
                                                                     bwa mem for HG00151
bwa mem for HG00106
bwa mem for HG01914
output:
                                                                      bwa mem for HG01985
bwa mem for HG01986
                                                                      bwa mem for HG02013
bwa mem for HG02051
                                                                      bwa mem for HG01879
bwa mem for HG01880
                                                                      bwa mem for HG01896
                                                                    bwa mem for NA18627 ]
                                                                                                             0
                                                                     bwa mem for HG01915
                                                                    bwa mem for NA18591
bwa mem for NA18566
                                                                                                             0
                                                                    bwa mem for NA18644
bwa mem for NA18545
                                                                                                             0 0 0
                                                                    bwa mem for HG00113
                                                                    bwa mem for HG00243
```

Now recursively search the present directory and all subdirectories for files with the .sam extension:

```
find . -name \*sam # the \* is a wild card to match anything in a file name before "sam"
```

```
/NA18591/NA18591.sam
/HG00260/HG00260.sam
/NA18909/NA18909.sam
/HG01985/HG01985.sam
/HG01879/HG01879.sam
/HG01915/HG01915.sam
/NA18627/NA18627.sam
/HG00106/HG00106.sam
/NA18644/NA18644.sam
 /NA19093/NA19093.sam
 /HG01880/HG01880.sam
/HG01986/HG01986.sam
,
/NA18907/NA18907.sam
/HG02013/HG02013.sam
/HG00149/HG00149.sam
/NA18757/NA18757.sam
/NA19098/NA19098.sam
/NA18870/NA18870.sam
/HG02051/HG02051.sam
/HG00243/HG00243.sam
/NA19137/NA19137.sam
/HG00151/HG00151.sam
/HG00132/HG00132.sam
/NA18545/NA18545.sam
/HG00113/HG00113.sam
/HG01896/HG01896.sam
/NA19256/NA19256.sam
/NA19138/NA19138.sam
 /HG01914/HG01914.sam
/NA18566/NA18566.sam
  w2470@gr062 task2]$
```

Q2.2. Please report your grep command and find commands and there outputs in your report. How many .sam files were produced? What do the exit statuses of the 30 subjobs indicate? [1 point].

answer: there are 30 .sam files were created, and the exit statuses are 0s.

##Task 3: SAM/BAM format and Samtools In this task, you will use Samtools to extract information from and manipulate a SAM/BAM file. You may wish to first watch the pre-recorded lecture on the SAM/BAM format as this provides the necessary background and review the following online documentation for Samtools and the SAM format.

-Samtools documentation: http://www.htslib.org/doc/samtools.html (http://www.htslib.org/doc/samtools.html) -SAM format specification: https://samtools.github.io/hts-specs/SAMv1.pdf (https://samtools.github.io/hts-specs/SAMv1.pdf)

Log into Greene and check out a compute node:

```
srun --time=4:00:00 --mem=4GB --pty /bin/bash
```

For questions in Task 3, you will conduct a series of operations on a BAM from the 1000 genomes project with paired-end reads mapped to chromosome 20. The path to the BAM file on Greene is:

```
/scratch/work/courses/BI7653/hw3.2022/HG00096.chrom20.ILLUMINA.bwa.GBR.low_coverage.20120522.bam
```

For all answers for Task 3, please report your command lines and answer the questions. You can find assistance in the above mentioned samtools and SAM specification links and in the course lecture(s) for this week. Q3.1. Review the samtools view documentation. Then, use this program to extract only the header from the bam file above and answer the following [1 point] Q3.1a. Report your command line

samtools view -H /scratch/work/courses/BI7653/hw3.2022/HG00096.chrom20.ILLUMINA.bwa.GBR.low\_coverage.20120522.bam

Q3.1b. Report the @HD header tag line. What does the information in this line indicate?

```
@HD VN:1.0 SO:coordinate
```

Answer: it means this bam file is sorted.

Q3.2. Use samtools view to answer the following. Review samtools view options -c, -f, and -F. Please answer the following questions including (1) your command line you used to obtain the answer and (2) the output written to your terminal [1 point].

Q3.2a how many unmapped reads are there in the BAM (hint; use appropriate bitwise flag(s) described in SAM/BAM lecture and documentation?

 $samtools\ view\ -c\ -f\ 4\ /scratch/work/courses/B17653/hw3.2022/HG00096.chrom20.ILLUMINA.bwa.GBR.low\_coverage.20120522.bam$  ## output: 7247

Q3.2b How many mapped reads are there in the BAM?

samtools view -c -F 4 /scratch/work/courses/B17653/hw3.2022/HG00096.chrom20.ILLUMINA.bwa.GBR.low\_coverage.20120522.bam ## output: 2924253

Q3.2c What is the percentage mapping rate (total mapped reads / total reads in the alignment) for this sample?

```
## 99.75% mapping rate
```

Q3.3. A hypothetical SAM file has alignment records with the bitwise flag values that include 4, 147, 113, 99 on the decimal scale. What are the binary and hexadecimal representations of each of the these values? [1 point]. Hint: https://www.rapidtables.com/convert/number/binary-to-hex.html)

```
Answer:
## Binary Hexadecimal
## 4 - 100 4
## 147 - 10010011 93
## 113 - 1110001 71
## 99 - 1100011 63
```

Q3.4. Picard has an online tool for determining the meanings of bitwise flag values such as those in Q3.3: http://broadinstitute.github.io/picard/explain-flags.html (http://broadinstitute.github.io/picard/explain-flags.html) Using this tool, what are the characteristics of the reads in each these four flags in Q3.3 (4.147,113,99)? [1 point].

```
Answer:
# # 4 - read unmapped
# # 147 - read paired, read mapped in proper pair, read reverse strand, scond in pair
# # 113 - read paired, read reverse strand, mate reverse strand, first in pair
# # 99 - read paired, read mapped in proper pair, mate reverse strand, first in pair
```

## ngs\_week3 assignment - Xiaoyan Wen

Q3.5. SAM specification allows for three types of alignment records. There are primary alignments, secondary alignments and supplementary alignments. Depending on the alignment software and command line used, secondary or supplementary alignments may also exist in a SAM/BAM file.BWA-MEM will typically add both primary and secondary alignments. Note that by using the -M option above we instructed bwa to set all supplementary alignments to secondary (0x100).

For Q3.5, use samtools view with appropriate -c, -f, -F options to count the following in the BAM file at

/scratch/work/courses/BI7653/hw3.2022/HG00096.chrom20.ILLUMINA.bwa.GBR.low\_coverage.20120522.bam

For each answer, provide the number of reads and the command line you used. Now answer the following [1 point]: Q3.5a How many alignments are primary?

samtools view -c -F 256 /scratch/work/courses/BI7653/hw3.2022/HG00096.chrom20.ILLUMINA.bwa.GBR.low\_coverage.20120522.bam ## output: 2931500

Q3.5b How many alignments are secondary?

samtools view -c -f 256 /scratch/work/courses/BI7653/hw3.2022/HG00096.chrom20.ILLUMINA.bwa.GBR.low\_coverage.20120522.bam ## output:  $\theta$ 

Q3.5c How many alignments are supplementary

samtools view -c -f 2046 /scratch/work/courses/BI7653/hw3.2022/HG00096.chrom20.ILLUMINA.bwa.GBR.low\_coverage.20120522.bam ## output: 0

Q3.5d What is the number of reads excluding unmapped reads, supplementary reads, secondary reads and PCR duplicates?

samtools view -c -F 3332 /scratch/work/courses/B17653/hw3.2022/HG00096.chrom20.ILLUMINA.bwa.GBR.low\_coverage.20120522.bam ## output: 2885340

Q3.6. A common task is to subset a SAM/BAM to include a subset of positions on a chromosome. Use samtools view to subset the BAM from Q3.5 from chromosome 20 position 1 to 2000000 (i.e., 2 Mb), while also retaining the header. Note that to perform this type of operation, the BAM must be coordinate-sorted, (which it is). Hint: See http://www.htslib.org/doc/samtools.html (http://www.htslib.org/doc/samtools.html)

Now count the reads in the subsetted BAM. Report the command line used to subset the BAM and the number of reads in the subset [ 1 point ].

samtools view -c /scratch/work/courses/BI7653/hw3.2022/HG00096.chrom20.ILLUMINA.bwa.GBR.low\_coverage.20120522.bam 20:1-20000

## output: 95338